by **eyeV**

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### Skyle (Remote Eyetracker Device)

---

Supply voltage	5V DC
Current consumption	0.9 A
Sampling rate	18 Hz (max. 30Hz)
Accuracy	1° - 2° (better accuracy planned)
Eye tracking mode	Binocular
Operating distance	45 - 65cm
Tracking Range (Width x Height)	30 x 12cm
Illuminators	NIR 850nm
Dimensions	200 x 42,5 x 44mm (l*w*h)

### System

---

Device Compatibility	iPad Pro 12.9 3 <sup>rd</sup> Generation (2018)
OS Compatibility	iOS 13

#### **Please note:**

Skyle for iPad uses the Assistive Touch function of the iPad. Before you update your iPad, make sure that the new version is compatible with eyetracker. For more information, please visit [eyeV.de](https://eyeV.de) or contact us at [support@eyeV.de](mailto:support@eyeV.de)

*eyeV GmbH reserves the right to make modifications, improvements, corrections or other changes without further notice to this document and any product described herein.*

*eyeV GmbH products are specifically not authorized for use as critical components in life support devices or applications, where the use can be reasonably expected to result in significant injury to the user.*

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*Settings -> Accessibility -> Touch -> AssistiveTouch -> Customise Top-Level Menu*

**Published by eyeV GmbH**  
**Alzeyer-Straße. 39, D-67549 Worms**  
**www.eyev.de © 2019 eyeV GmbH. All rights reserved.**

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

*Settings -> Accessibility -> Touch -> AssistiveTouch -> Dwell Control*

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*Settings -> Accessibility -> Touch -> AssistiveTouch -> Movement Tolerance*

**Published by eyeV GmbH**  
**Alzeyer-Straße. 39, D-67549 Worms**  
**www.eyev.de © 2019 eyeV GmbH. All rights reserved.**

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
*Settings -> Accessibility -> Touch -> AssistiveTouch -> Custom Gestures -> Create New Gesture...*

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
*Settings -> Home Screen & Dock -> App Icons*

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*Settings -> Display & Brightness -> Display Zoom -> View -> Zoomed*

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